

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Decellular Dermis Cryo
Product code(s): 04 - 351
Batch Number(s): 83210301
Date Manufactured: 15 Mar 2021
Expiry Date: 14 Mar 2022
Appearance: Clear Solution

Quality Control:

	Specification	This batch result
pH	7.10 - 7.50	7.23
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to 430/3 Effective 19OCT2017

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter and filled into Gamma Irradiated Bottle

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *S.J. Jaylor*

Date: 22 MAR 21

Signed for Certificate Check: *T Waller*

Date: 22 MAR 21