

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Low Ionic Strength Saline (LISS)

**Product code(s):** 03210301

**Batch Number(s):** 04-016

**Date Manufactured:** 29 Mar 21

**Expiry Date:** 28 Mar 22

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.50 - 6.90	6.82
Conductivity( $\mu\text{s}/\text{cm}$ )	3500 - 3800	3643
Osmolality(mOsm/kg)	285 - 295	292

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** Sodium azide (0.01% [w/v])

**Sterility:** N/A

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 30 MAR 21

**Signed for Certificate Check:**

**Date:** 30 mar 21