

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Frozen Cell Recovery Solution (FCRS)

Product code(s): 04 - 161 BAG

Batch Number(s): 41210301

Date Manufactured: 22 Mar 2021

Expiry Date: 21 Mar 2023

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.40 - 6.80	6.58
Conductivity($\mu\text{s}/\text{cm}$)	44000 - 52000	48240
Osmolality(mOsm/kg)	1080-1120	1108

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: EDTA (anti-fungal)

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *S. J. Taylor*

Date: 29 MAR 21

Signed for Certificate Check: *T. Waller*

Date: 29 MAR 21