

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** AmnioCryo Medium (Amniotic Membrane Cryopreservation)  
**Product code(s):** 04-326  
**Batch Number(s):** 72210401  
**Date Manufactured:** 08 Apr 2021  
**Expiry Date:** 07 Apr 2023  
**Appearance:** Clear Red Solution  
**Quality Control:**

|    | Specification | This batch result |
|----|---------------|-------------------|
| pH | 7.10-7.50     | 7.44              |

Manufactured to SPN431/3 – Effective date 21 July 2020.

**Preservatives used:** N/a

**Sterility:** Product is filtered through a 0.2µm filter and filled into gamma irradiated bottles..

**Storage:** Ambient.

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 14 APR 21

**Signed for Certificate Check:**



**Date:** 14 APR 21