

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Cambridge Antibiotic Solution

**Product code(s):** 04-301

**Batch Number(s):** 70210401

**Date Manufactured:** 12 Apr 2021

**Expiry Date:** 31 Jan 2022

**Appearance:** Red Suspension

**Quality Control:**

	Specification	This batch result
pH	7.1 - 7.5	7.30
Osmolality(mOsm/kg)	280 - 300	297

Manufactured to SPN 310/3 effective 07 May 19

**Preservatives used:** Gentamycin Sulphate, Imipenem, Polymixin B, Vancomycin, Nystatin

**Sterility:** Product is filtered through a 0.2µm filter prior to addition of Nystatin in a clean room and filled into Gamma Irradiated bottles.

**Storage:** <-40 °C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 16 April 2021

**Signed for Certificate Check:** 

**Date:** 16 APR 2021