

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsevers Solution

**Product code(s):** 04 - 182HE

**Batch Number(s):** 32210401

**Date Manufactured:** 20 Apr 2021

**Expiry Date:** 19 Apr 2022

**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.04
Conductivity( $\mu\text{s}/\text{cm}$ )	11500 - 13000	12170

Manufactured to SPN 870/2.1 Dated 18.08.16

**Preservatives used:** Chloramphenicol and Neomycin Sulphate

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2 - 8°C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 28 Apr 21

Signed for Certificate Check: 

Date: 28 APR 21