

CERTIFICATE OF ANALYSIS

## Certificate of Analysis

**Product name:** Modified Alsevers Solution

**Product code(s):** 04 - 182HE

**Batch Number(s):** 32210402

**Date Manufactured:** 26 Apr 2021

**Expiry Date:** 25 Apr 2022

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.80 - 7.30	7.06
Conductivity( $\mu\text{s}/\text{cm}$ )	11500 - 13000	12250

Manufactured to SPN 870/2.1 Dated 18.08.16

**Preservatives used:** Chloramphenicol and Neomycin Sulphate

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2 - 8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 30 Apr 21

**Signed for Certificate Check:** 

**Date:** 30 Apr 21