

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution

Product code(s): 04 - 182HE

Batch Number(s): 32210601

Date Manufactured: 14 June 2021

Expiry Date: 13 June 2022

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.05
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	12147

Manufactured to SPN 870/2.1 Dated 18.08.16

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8°C


Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 21 June 21

Signed for Certificate Check:



Date: 21 Jun 21